
Product data sheet





<u>+44 (0)1946 816600</u>



SPH011 HYDROMETER - SPECIAL PURPOSE HYDROMETER - LATEXOMETER

These instruments are calibrated in grams of dry rubber per litre of latex and adjusted at 28.9°C (84°F). 250mm long.

All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product features:

- · Calibrated in grams of dry rubber per litre of latex
- Adjusted at 28.9°C (84°F)
- 250mm long
- · Made in UK
- · Supplied in protective, resealable plastic tube

Product specifications:

Material: Soda glass, steel shot ballast, synthetic wax resin

Accuracy: +/- 10.0

Ranges: 50 to 450 g/l or 50 to 250 g/l

26mm or 43mm Bulb diameter:

Overall length: 250mm

Standard product numbers:

Range	Bulb diameter	Gross product weight*	Product number
50 to 450 g/l	26mm	100g	10/611/0
50 to 250 g/l	43mm	133g	10/612/0

^{*}including resealable plastic tube



View me on the website



